

Bionexis

Laboratory

Phone 9310733411

Email info@bionexis.in

Web. www.bionexis.in

Address:- Plot No.A-28, Mohan Co-operative Industrial Estate,
Mathura road, New Delhi, New Delhi, Delhi, 110044



BND0001
DO Meter

Type	MICROPROCESSOR BASED
Function	DISSOLVED OXYGEN / TEMPERATURE SIMULTANEOUSLY DISPLAY
DO Range	0 to 20 ppm
Temp Range	0 to 70° C
Resolution	DO 0.01 ppm temp 0.1° C
Accuracy	±1% ppm
Display	16x2 Alphanumeric
Measuring Cell	Amperometric Gold / Silver Electrode
Temperature Compensation	Automatic
Memory	Storage of 100 sample
Dimensions, Weight	275 x 187 x 94mm (approx), 1.5 kgs (approx)
Accessories	DO cell cell stand temp. Probe Operating manual
Price	9,322/- Prices are inclusive of all taxes, Shipping and Packaging charges 1 Year Warranty*, 10 Days Replacement, COD is Available