

Bionexis

Laboratory

Phone 9310733411

Email info@bionexis.in

Web. www.bionexis.in

Address:- Plot No.A-28, Mohan Co-operative Industrial Estate,
Mathura road, New Delhi, New Delhi, Delhi, 110044



BNPE001

PH Electrode ,acrylic unbreakable Body Ph Range(0 to 14)

| | |
|-------------------------|---|
| PH Range | 0 to 14 |
| Liquid Junction | Ceramic |
| Cable | 1m |
| Connector | BNC |
| Overall Length | 18 cm |
| Filling Solution | 3.33M KCl |
| Body Material | acrylic |
| Tip / Shape | spheric |
| Body | Acrylic |
| Price | 826/- Prices are inclusive of all taxes, Shipping and Packaging charges 10 Days Replacement COD is Available |