

Bionexis

Laboratory

Phone 9310733411

Email info@bionexis.in

Web. www.bionexis.in


Address:- Plot No.A-28, Mohan Co-operative Industrial Estate,
Mathura road, New Delhi, New Delhi, Delhi, 110044



Smart Ph Meter

PH Range	0.00 ph to 14.00 ph
Display	LCD, 16x2 Alphanumeric
Range	0.00 ph to 14.00 ph , 0 to \pm 1999.9mv , temp -20.0 to 150° C
Resolution	0.01 pH , 0.1mv , 0.1° C
Accuracy	\pm 0.05%
Temperature Compensation	Automatic/ Manual
Auto Buffer Recognition	4.00 , 7.00 , 9.20
Calibration	3 point
Dimensions	275 x 187 x 94mm (approx)
Power	220V \pm 10V AC
Accessories	<ol style="list-style-type: none"> 1. pH meter, 2. pH Electrode 3. Temperature Sensor (Not available in BNPH001) 4. Stand 5. 3 Buffers : 4.00 , 7.00 , 9.20 6. Operating manual

What's In the Box



pH Meter



Stand



pH Electrode



Standard Buffers



Temperature Sensor



User Manual

#Cat No	Description	Price
BNPH001	Smart Ph Meter, Manual Temperature Compression	5,605.00
BNPH002	Smart Ph Meter With Automatic Temperature Compression	7,021.00
BNPH003	Smart Ph Meter With Automatic Temperature Compression and Battery Backup	8,791.00

Prices are inclusive of all taxes, Shipping and Packaging
 1 Year Warranty*
 10 Days Replacement
 COD is Available