

Bionexis

Laboratory

Phone 9310733411

Email info@bionexis.in

Web. www.bionexis.in

Address:- Plot No.A-28, Mohan Co-operative Industrial Estate,
Mathura road, New Delhi, New Delhi, Delhi, 110044



Microprocessor Based TDS/Conductivity/Temp. Meter

Function	Display TDS / EC / TEMPERATURE SIMULTANEOUSLY
TDS Range	0 to 1000 ppt Auto ranging
Electronic Conductivity Range	0 to 1000 ms Auto ranging
Temperature	0 to 100° C
Resolution	TDS 0.1 ppm
Electronic Conductivity	0.1 μ mhos
Temperature	0.1 degree C
Accuracy	0.5% of range
Cell Constant	Adjustable
Display	LCD, 16x2 Alphanumeric
Measuring Cell	Platinum DIP Type
Temperature Compensation	Automatic
Dimensions	275 x 187 x 94 mm (Approx)
Weight	1kg (approx)
Power Supply	230 V \pm 10 V AC , 50Hz
Accessories	Conductivity cell Stand Temperature Sensor (NA in BNTDS001) Operating manual

What's In the Box



#Cat No	Description	Price
BNTDS001	TDS/Conductivity meter without automatic temperature compression	6,903.00
BNTDS002	TDS/Conductivity/Temp. Meter with automatic temperature compression	7,375.00

Prices are inclusive of all taxes, Shipping and Packaging charges
 1 Year Warranty*
 10 Days Replacement
 COD is Available