

# Bionexis

## Laboratory

**Phone** 9310733411

**Email** [info@bionexis.in](mailto:info@bionexis.in)

**Web.** [www.bionexis.in](http://www.bionexis.in)

**Address:-** Plot No.A-28, Mohan Co-operative Industrial Estate,  
Mathura road, New Delhi, New Delhi, Delhi, 110044




BNPHS001

## Smart pH Meter with Inbuilt Magnetic Stirrer

<b>PH Range</b>	0.00 ph to 14.00 ph
<b>Display</b>	16x2 Alphanumeric
<b>Range</b>	0.00 ph to 14.00 ph , 0 to $\pm 1999.9\text{mv}$ , temp -20.0 to 150° C
<b>Resolution</b>	0.01 pH , 0.1mv , 0.1° C
<b>Accuracy</b>	$\pm 0.05\%$
<b>Temperature Compensation</b>	Automatic/ Manual
<b>Auto Buffer Recognition</b>	4.00 , 7.00 , 9.20
<b>Calibration</b>	3 point
<b>Dimensions</b>	275 x 187 x 94mm (approx)
<b>Power</b>	220V $\pm$ 10V AC
<b>Accessories</b>	<ol style="list-style-type: none"> <li>1. pH meter with magnetic stirrer,</li> <li>2. pH Electrode</li> <li>3. Temperature Sensor</li> <li>4. Stand</li> <li>5. 3 Buffers : 4.00 , 7.00 , 9.20</li> <li>6. 100ml Beaker</li> <li>7. Operating manual</li> </ol>

# What's In the Box



**pH Meter with  
Magnetic Stirrer**



**pH Electrode**



**Temp. Sensor**



**Stand**



**Standard Buffers**



**100ml Beaker**



**Magnetic Bead**



**User Manual**

#Cat No	Description	Price
<b>BNPHS001</b>	Ph Meter With Automatic Temperature Compensation and inbuilt Magnetic Stirrer	9,086.00
Prices are inclusive of all taxes, Shipping and Packaging 1 Year Warranty* 10 Days Replacement COD is Available		