

Bionexis

Analytical Instruments

Phone 9310733411

Email info@bionexis.in

Web. www.bionexis.in

Address:- Plot No.A-28, Mohan Co-operative Industrial Estate,
Mathura road, New Delhi, New Delhi, Delhi, 110044

Index

Analytical Instrument

1. Ph Meters
2. Multiparameter
3. TDS/Conductivity Meters
4. Dry Bath Incubators
5. Haemoglobin Meter
6. Turbidity Meter
7. Photocolorimeter
8. DO Meter
9. Accessories



Smart Ph Meter

PH Range	0.00 ph to 14.00 ph
Display	LCD, 16x2 Alphanumeric
Range	0.00 ph to 14.00 ph , 0 to ± 1999.9mv , temp -20.0 to 150° C
Resolution	0.01 pH , 0.1mv , 0.1° C
Accuracy	± 0.05%
Temperature Compensation	Automatic
Auto Buffer Recognition	4.00 , 7.00 , 9.20
Calibration	3 point
Dimensions	275 x 187 x 94mm (approx)
Power	220V ± 10V AC
Accessories	1. pH meter, 2. pH Electrode 3. Temperature Sensor (Not available in BNPH001) 4. Stand 5. 3 Buffers : 4.00 , 7.00 , 9.20 6. Operating manual

What's In the Box



pH Meter



Stand



pH Electrode



Standard Buffers



Temperature Sensor



User Manual

#Cat No	Description	Price
BNPH002	Smart Ph Meter With Automatic Temperature Compression	6842.37
BNPH003	Smart Ph Meter With Automatic Temperature Compression and Battery Backup	8567.50
• 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement		



BNPHS001
Smart pH Meter with Inbuilt Magnetic Stirrer

PH Range	0.00 ph to 14.00 ph
Display	16x2 Alphanumeric
Range	0.00 ph to 14.00 ph , 0 to ± 1999.9mv , temp -20.0 to 150° C
Resolution	0.01 pH , 0.1mv , 0.1° C
Accuracy	± 0.05%
Temperature Compensation	Automatic
Auto Buffer Recognition	4.00 , 7.00 , 9.20
Calibration	3 point
Dimensions	275 x 187 x 94mm (approx)
Power	220V ± 10V AC
Accessories	1. pH meter with magnetic stirrer, 2. pH Electrode 3. Temperature Sensor 4. Stand 5. 3 Buffers : 4.00 , 7.00 , 9.20 6. 100ml Beaker 7. Operating manual

What's In the Box




#Cat No	Description	Price
BNPHS001	Ph Meter With Automatic Temperature Compensation and inbuilt Magnetic Stirrer	8855.00
• 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement		



BNMP001
 Laboratory Multiparameter Analyser pH/ORP/TDS
 /Conductivity/Temperature

PH Range	0.00 ph to 14.00 ph
PH Resolution	0.02
PH Accuracy	±0.02
Calibration Points	Three
mv(ORP) Range	-500.0 to 500.0 mv
mv(ORP) Resolution	0.2 mv
TDS Range	0 to 5 ppt Auto Ranging
TDS Resolution	0.1 ppm
TDS Accuracy	1% of Range
Conductivity Range	0 to 10 ms
Conductivity Resolution	0.1 us
Calibration Points	One
TDS Results	Conductivity Conversion
Cell Constant	1.00
Temperature Range	5 to 70° C
Temperature Compensation	Automatic

What's In the Box



Display	16x2 LCD
Power	220 Volt $\pm$ 10%, 50 Hz
Dimensions	275x187x94 mm (Approx)
Accessory	1. Multiparameter Analyser 2. Electrode Stand, 3. pH Electrode, 4. Conductivity Cell 5. Buffer Tablets: 4.00, 7.00 and 9.20, 6. Instruction Manual

#Cat No	Description	Price
BNMP001	Multiparameter pH/ORP/TDS/Conductivity/Temperature	10925.00
• 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement		



Microprocessor Based TDS/Conductivity/Temp. Meter

Function	Display TDS / EC / TEMPERATURE SIMULTANEOUSLY
TDS Range	0 to 1000 ppt Auto ranging
Electronic Conductivity Range	0 to 1000 ms Auto ranging
Temperature	0 to 100° C
Resolution	TDS 0.1 ppm
Electronic Conductivity	0.1 μ mhos
Temperature	0.1 degree C
Accuracy	0.5% of range
Cell Constant	Adjustable
Display	LCD, 16x2 Alphanumeric
Measuring Cell	Platinum DIP Type
Temperature Compensation	Automatic
Dimensions	275 x 187 x 94 mm (Approx)
Weight	1kg (approx)
Power Supply	230 V $\pm$ 10 V AC , 50Hz
Accessories	Conductivity cell Stand Temperature Sensor (NA in BNTDS001) Operating manual

What's In the Box



TDS Meter



Conductivity Cell



Temp. Sensor



Stand



User Manual

#Cat No	Description	Price
BNTDS002	TDS/Conductivity/Temp. Meter with automatic temperature compression	7187.5
• 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement		



BNDI001

Dry Bath Incubator with 24 Test Tube Block

Display	16x2 Alpha Numeric LCD
Temperature Range	Ambient to 70° C
Temperature Resolution	0.1 C
Inbuilt Timer	0 to 99 Minute Timer
Power Supply	220 V $\pm$ 10V AC
Weight	1.5 kg Approx.
Dimensions	275 x 187 x94mm (Approx)
Accessories	Instruction Manual

#Cat No	Description	Price
BNDI001	Dry Bath Incubator 24 Block	8107.5
• 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement		



BNHB001
 Microprocessor Based Smart Haemoglobin Meter,
 Cyanmethoglobin Method

Measurement Method	Cyanmethoglobin Method
Display	16 x 2 Alphanumeric display
Wavelength	540nm
Calibration	automatic by press CAL key
Set Zero	automatic by push switch
Measuring Range	0.0 to 30.0 g / dl
Photo Detector	silicon photodiode
Accuracy	±4 %
Light Source	540nm green LED
Power	220 volt ± 10 AC
Frequency	50Hz
Dimensions	275 x 187 x94mm (approx)
Cuvette	12x 65 mm round cuvette
Weight	0.8 kg (approx)
Accessories	5 Cuvette, 1 Cuvette Holder, Operation manual
Set Of Round Cuvettes	5 Nos.
Test Tube Holder	01 No.
Price	8440.00



BNTB001

Microprocessor Based Turbidity Meter

Type	MICROPROCESSOR BASED
Range	0 to 1000 NTU
Resolution	1 NTU
Sample System	±2% FS
Light Source	high power white LED
Detector	photodiode
Display	16x2 Alphanumeric
Measuring Cell	platinum DIP type
Temperature Compensation	automatic
Dimension	275 x 187 x 94 mm (approx)
Weight	1 kg (approx)
Power Supply	220 V ± 10V AC , 50 Hz
Accessories	Test tube set Operating manual
Price	7072.00 • 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement



BNPC001

Photometer

Set Zero (Auto Zero)	by press push switch with audio response
Wavelength Range	405 - 700 nm
Filters	built - in 8 digital filters > 30,000 hours life
Filters Wavelength	405nm , 450nm , 480nm , 520nm 540nm , 578nm , 620nm , 700nm
Display Parameters	O.D, % T, wavelength
Detector	silicon photo diode
Accuracy	+ 0.01 O.D.
Stability	+0.02 O.D per hour
Light Source	highly accurate LED source
Volume	1ml minimum
Dimension	275 x 187 x94mm (approx)
Power	230 $\pm$ 10 V , 50 Hz
Accessories	set of test tube - 5 nos. test tube holder 1 no. Operating manual
Price	5461.00 Prices are inclusive of all taxes, Shipping and Packaging charges 1 Year Warranty* 10 Days Replacement COD is Available



BNDO001
DO Meter

Type	MICROPROCESSOR BASED
Function	DISSOLVED OXYGEN / TEMPERATURE SIMULTANEOUSLY DISPLAY
DO Range	0 to 20 ppm
Temp Range	0 to 70° C
Resolution	DO 0.01 ppm temp 0.1° C
Accuracy	±1% ppm
Display	16x2 Alphanumeric
Measuring Cell	Amperometric Gold / Silver Electrode
Temperature Compensation	Automatic
Memory	Storage of 100 sample
Dimensions, Weight	275 x 187 x 94mm (approx), 1.5 kgs (approx)
Accessories	DO cell temp. Probe Operating manual
Price	9084.00 • 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement



BNPE001
 PH Electrode ,acrylic unbreakable Body Ph Range(0 to 14)

PH Range	0 to 14
Liquid Junction	Ceramic
Cable	1m
Connector	BNC
Overall Length	18 cm
Filling Solution	3.33M KCl
Body Material	acrylic
Tip / Shape	spheric
Body	Acrylic
Price	805.00 • 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement



BNTS001
Temperature Probe

Measurement range	0 to 100° C
Accuracy	+/- 0.8C
Cable	75 cm
Connector	Pin
Overall Length	19 cm
Filling Solution	3.33M KCl
Body Material	Stainless Steel
Warranty	1 Year
Price	747.00 • 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement



BNST001
Stand Model 1

Price	• 559.00 • 18% GST Extra • Free Shipping • 1 Year Warranty • 10 Days Replacement
-------	--



BNST002
Stand Model 2

Price	919.00
-------	--------

Terms & Conditions

Shipping Policy

- Orders are processed within 2 business days and shipped via our trusted logistics partners.
- Estimated delivery time varies based on location:
 - Domestic: 4-7 days
 - International: 7-12 days
- Customers will receive a tracking number once the order is shipped.
- Shipping charges are free on selected products.

Warranty Registration

- After purchasing a product, customers must register their warranty by sending the machine name and serial number via:
 - **WhatsApp:** 9310733411
 - **Email:** info@bionexis.in

Warranty & Aftersales Service

- We offer a **12-month warranty** on our products against manufacturing defects.
- Warranty does not cover damages due to improper handling, modifications, or unauthorized repairs.
- In case of any service requirements, customers must send the machine to our **Bionexis Service Center** at:

Nexis International

25/369 B Block, New Motinagar, Baba Karam Singh Marg, Delhi - 110015

Phone: 9310733411

Email: info@bionexis.in

For further assistance, reach out to our customer support at **9310733411** or **info@bionexis.in**.